

Sequencing Quality Control Report

Genomics Unit

March 26, 2014

1 Flowcell Information

Flowcell BC3UG0ACXX
Run Start 2014-03-14
Read Length 100 paired end
HiSeq Sequencing Kit Truseq v3


Sequencer Software HiSeq Control Software 2.0
Gerald version 1.8

2 Samples Information


LabID	Name	Index	NumberOfReadPairs
6491	CTRL	ATCACG	135218638
6492	RED	ACTTGA	144586011
6493	SMU1	TAGCTT	129219623
6494	MFAP	GGCTAC	131737035

3 Quality Plots


Sample 6491: CTRL



Sample 6492: RED



Sample 6493: SMU1



Sample 6494: MFAP

