

Sequencing Quality Control Report

Genomics Unit

June 26, 2014

1 Flowcell Information



Flowcell BC4WCWACXX
Run Start 2014-06-19
Read Length 50 single end
HiSeq Sequencing Kit Truseq v3

Sequencer Software HiSeq Control Software 2.0
Gerald version 1.8

2 Samples Information

LabID	Name	Index	NumberOfReads
7302	B cells K27ac	TTAGGC	42958853
7303	B cells K9me1	GATCAG	53147581
7304	input	ACTGAT	57277961

3 Quality Plots



Sample 7304: input

