

Sequencing Quality Control Report

Genomics Unit

December 25, 2015

1 Flowcell Information

Flowcell AC88UTACXX
Run Start 2015-12-23
Read Length 50 single end
HiSeq Sequencing Kit TruSeq v3


Sequencer Software HiSeq Control Software 2.2.58
Gerald version 1.8

2 Samples Information


LabID	Name	Index	NumberOfReads
12372	Klf 4	CTTGTA	63799107
12372	Klf 4	GCCAAT	60360610
12372	Klf 4	TTAGGC	45804028

3 Quality Plots

Sample 12372_CTTGTA:



Sample 12372_GCCAAT:



Sample 12372_TTAGGC:

95.5 % of base calls \geq Q30

